

# Prevent Falls & Wandering

## Early-Warning Chair Belts Know before resident gets up!



### Works with multiple chair exit alarms

#### Early Warning

The Fall Prevention Wheelchair Seat Belts reduce falls by triggering a fall monitor when the Easy Release Buckle is unbuckled or the Quick Release Hook & Loop strap is opened; notifying caregivers before a resident attempts to leave their chair. These early warning seat belts are designed for use with all of our TL-2100 Series Fall Monitors.

#### Non-Restraining

Both the Easy Release & Quick Release seat belts preserve resident's freedom of movement by not acting as a restraint. Further, the Quick Release Seat Belt is ideal for residents with low dexterity as the Hook & Loop strap can be easily unfastened with one hand.

**IMPORTANT:** Seat Belts are not designed to restrain or hold individuals in position in their chairs.



Easy Release TL-2109



Quick Release TL-2109V

#### Key Features:

- Easy Release Seat Belt features a large, red button that is easy to locate and release with one finger
- Quick Release Seat Belt features an easily grasped Hook & Loop strap which unfastens quickly even for those with limited dexterity
- **Fully Adjustable Belt Length** - 20" to 54"
- Metal ends install securely on most wheelchair models



**Call Us Today! (800) 650-3637 | [www.smartcaregiver.com](http://www.smartcaregiver.com)**